

Burton Outpatient® II Lights | OP216SC



Stock Photo – Options May Vary

FACTORY NEW

Burton OP216SC Outpatient® II Lights. The Outpatient® is the documented industry leader for procedure lights. More physicians choose the Outpatient® II than any other procedure light for focusable flood-type light, maximum intensity, and shadow-free operation, with minimal heat output and high color temperature.

Perfect for the hospital, surgical suite, or physician's office. Outpatient® II Single Ceiling Light Features: Focuses with central, removable, autoclavable SteriHandle™. The friction knob allows the user to adjust friction on the self-balancing arm positioning to the desired tension without a tool. Additionally, cool operation is assured with heat filters and dichroic coatings.

FEATURES

- Outpatient® II is the documented industry leader in procedure lights
- Focusable flood-type light, maximum intensity, and shadow-free operation
- Minimal heat output and high color temperature
- Perfect low-maintenance, top-performing light for nearly every application in the hospital, surgical suite, or physician's office
- Focuses with central removable, autoclavable SteriHandle™

SPECIFICATIONS

- 80,000 lux at 24" (61 cm)
- 3300 K color temperature
- CRI (Color Rendering Index): 94
- Removable, autoclavable SteriHandle™
- Friction knob

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)