

Gendron 890 Stretcher – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Our Completely Refurbished Gendron 890 Stretcher has all the features you need at a great value! The Easy Access Gendron 890 Stretcher features Trendelenburg and Reverse Trendelenburg, Pneumatic Assist Head Up/Down, and Easy Height Up/Down, with conveniently located foot pedals on either side.

Additional Gendron 890 Features: Central steering/braking system can be operated from either end/side. The Fowler-assisted Backrest raises 90° – when in the uprights position, conveniently allows access to the storage area, which houses an IV Pole (sold separately), security cable, and an O2 Tank Well.

Full-length Siderails fold down and under out of the way to allow easy access.

FEATURES

- Trendelenburg/Rev Trendelenburg
- Head-Up/Down (Pneumatic Assist)
- Height Up/Down
- Fold Under Siderails
- 4-Wheel Locking

SPECIFICATIONS

- 4 Wheel Locking
- 28"x76" patient surface
- 6 Months Parts Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)