

Amsco M70 Blanket Warmer / Warming Cabinet - Refurbished | M70WC-EL



Stock Photo - Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Amsco M70 is a cost-effective solution for blankets and fluid warming needs in your medical facility. The unit is Completely Refurbished to Excellent Operational and Cosmetic Condition. Analog controls for the top compartment to adjust the temperature. The bottom cabinet stays at a constant of 120°. Additional Amsco M70 Features: This Unit has Audible as well as visual alarms in case of over-temp detection.

Comes with a hospital-grade power cord (new). Also included a blanket rack in the top compartment, a shelf in the lower compartment, and clips. See below for more information.

Note: The Amsco M70 Blanket Warmer has been discontinued. Call for availability.

FEATURES

- Dual Compartment
- Free Standing Unit
- Top and Bottom Compartment Heated
- A shelf in the Bottom Compartment
- Stainless Steel Construction

SPECIFICATIONS

- Construction: Stainless Steel
- Voltage: 110v Completely Refurbished and Biomedically Inspected

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com